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009151EX-20150401

Oxyd-Stop macro für NEM + Modellguss

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Instructions for use

Oxyd-Stop macro for non-precious and chrome cobalt alloys

Before using the product, please read these instructions for careful use.

The manufacturer will not accept any liability for damage resulting from non-compliance with these instructions for use.

1. Description of the product

Oxyd-Stop macro is a ready-to-use solution with an integrated brush for the application to non-precious and CoCr alloys. A fine ceramic layer is obtained which avoids oxidation of most CoCr alloys.

2. Indication

Oxyd-Stop macro is applied to the polished or oxide-free surface and effectively avoids oxidation during heat treatment. Additionally, Oxyd-Stop macro can also be used as an anti-flux. It does not contain any flux; therefore it can be used for any type of ceramic restorations in all furnaces.

3. Contraindication

Based on the current state of knowledge there are no known contraindications.

4. Hazard and safety information

Signal word:

Danger

Hazard statements:

H225	Highly flammable liquid and vapour.
H319	Causes serious eye irritation.
H336	May cause drowsiness or dizziness.
H361d	Suspected of damaging the unborn child.

Precautionary statements:

P201	Obtain special instructions before use.
P210	Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.
P233	Keep container tightly closed.
P261	Avoid breathing vapours.
P280	Wear protective gloves / eye protection.
P281	Use personal protective equipment as required.

Response:

P303+P361+P353	IF ON SKIN (or hair): Remove / Take off immediately all contaminated clothing. Rinse skin with water / shower.
P305+P351+P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P308+P313	IF exposed or concerned: Get medical advice / attention
-----------	---

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Storage:

P403+P235	Store in a well-ventilated place. Keep cool.
-----------	--

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Shake Oxyd-Stop macro well before use since it sediments after a short time. Apply a covering coat to the metal surface and allow to dry for approx. 10 seconds. Then the metal can be heated.

Oxyd-Stop macro should not be applied to the ceramic surface to avoid the formation of rough surfaces; however, a chemical reaction with the ceramic will not occur.

Sandblast residual Oxyd-Stop macro with glass beads at a maximum pressure of 3 bars.

If Oxyd-Stop macro is no longer liquid after it is left standing over a longer period or cannot be shaken, it can be diluted using the thinner, REF 520 0064 0.

7. Symbols

7.1 Symbols according to DIN EN ISO 15223-1

	Reference number
	Lot number
	Can be used until YYYY-MM
	Manufacturer

7.2 Symbols according to GHS regulation

	GHS02, flame
	GHS07, hazardous to health
	GHS08, systemic health hazards

7.3 Symbols according to NFPA (USA)

	Hazard diamond
--	----------------

8. Additional information

The product may only be used by dentists, dental technicians and adequately trained dental staff.

Only original tools and components may be used for processing. The lot numbers of all components should be recorded to ensure traceability and to handle complaints.

These instructions for use are based on state-of-the-art methods and equipment and our own experience. The product may only be used for the indication described under item 2. The user himself is responsible for processing the product. Liability for incorrect results shall be excluded since the manufacturer does not have any influence on processing. Any occurring claims for damages may only be made up to the value of our products.

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